



Colony Counter 220-240V 50Hz 1Ph.

Code: [535-210/01](#)

Product Group: [Microbiology](#)

The petri dish is placed on the viewing and every time a colony is pressed a count is registered and accompanied by an audible tone. Both the amount of pressure needed to register a count and the volume of the tone are adjustable. The built-in adjustable magnifier enables samples to be examined more easily. For use with up to 100 mm glass or plastic petri dishes.

- Uniform glare-free illumination.
- Angled front for operator comfort.
- Counter triggered by depressing petri dish.
- Digital display of results.
- Built in magnifier.

Further Information

For 220-240 V, 50 Hz, 1 ph.