



DS7.2 Direct and Residual Shear Strength Program for Windows 7, 32/64 bit

Code: 27-1793

Product Group: DS7, DS7.2 Direct/Residual Shear Software, Automatic Data Acquisition - Direct Shear

This unique package provides test options for quick undrained or drained shear tests with the user selectable option of residual testing. Individual test results can be linked together to produce the Coulomb Envelope. Printouts and plots are available for sample description and basic test data such as moisture content, etc. Realtime plots of settlement, shear versus displacement and vertical displacement during shearing is readily available via the PC screen or printer.