



Laboratory Vane Apparatus with 4 Springs and 12.7mm Vane; Hand Operated.

Code: [26-2270](#)

Product Group: [Laboratory Vane Apparatus](#),
[Laboratory Vane Apparatus](#)

The equipment is based on a design by the Transport Research Laboratory, England, and is available in hand or motorised versions. Stress is applied through the 12.7 x 12.7 mm vane by means of any one of four calibrated springs.

Standards

BS 1377

Further Information

Hand-operated, complete with 4 calibrated springs and a 12.7 mm diameter x 12.7 mm long vane.

Accessories




Attachment Clamp for 38 mm and 100 mm Sample Tubes

Code: [26-2289](#)



Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)

Spares/Consumables



Laboratory Vane Apparatus (Motorised)

Code: [26-2275/01](#)



Set of 4 Calibrated Springs for Laboratory Vane Apparatus.


Code: [26-2275/10](#)

Alternatives



Attachment Clamp for 38 mm and 100 mm Sample Tubes

Code: [26-2289](#)




Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)