



## Spore Sampler

Code: [535-260/01](#)

Product Group: [Microbiology](#)

- Continuous sampling for up to 7 days.
- Built-in vacuum pump.
- 7 jewel clockwork recording drum.

The spore sampler uses a pump to draw air through an orifice. Fungal spores and pollen are impacted onto a strip of adhesive coated plastic strip fitted to the clockwork driven drum. The sampling head can revolve and a large fin is fitted so that air is always drawn from upwind.

Complete with Sampler, Wind Vane, Rain Sheild, Drum, Cutting Block, 1 kg Gelvatol, Tape and Scalpel

## Further Information

220-240 V, 50 Hz, 1 ph